

Serological Pipette

DATA SHEET



ISO9001: 2015 ISO13485: 2016




Product Name: Serological Pipette	
Product Materials:	Pipette: Polystyrene (PS).
	Filter: Polyethylene (PE).
Product Shelf Life:	Temperature Range: Room temperature storage.
	Shelf Life: 3 years after date of production (provided the packaging remains intact).
Sterilization:	Yes.
	Sterilized by electron beam irradiation, SAL=10 ⁻⁶ .
	This product has been irradiated and released based on ANSI/ AAMI/ ISO 11137 standards.
Biological Safety & Purity	Free of detectable DNase/RNase
	Non-Pyrogenic
	Non-Cytotoxic
	BSE/TSE free
Performance Test:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures.
	Appearance inspection: qualified Stress test: qualified
	Sealing test: qualified Packaging inspection: qualified
Features:	1. Made from high clarity polystyrene.
	2. Withstand temperature range: Approx. 0°C to 40°C.
	3. Individual packaging.
	4. Built-in filter element effectively prevents potential cross-contamination.
	5. Clear and accurate bidirectional graduation lines.
	6. Negative scale design provides extra capacity.
	7. Color-coded rings for easy capacity identification.
	8. Sterilized by E-beam.
	9. DNase/RNase/Human DNA free.
	10. Non-pyrogenic, non-cytotoxic.
Product Application:	Serological pipettes are essential tools for accurately transferring measured volumes of liquid in cell culture, bacteriological, and research laboratories.

DATA SHEET

Technical Drawing


1 mL Standard Serological Pipette



Stretch

Item	Product Name	Graduation/mL	Volume	Sterile	Packaging	Packaging Specification
J2001	1 mL Standard Serological Pipette	1/100	1mL	Yes	Individually paper-plastic wrapped	50 Pieces/Bag, 20 Bags/Case
J2201	1 mL Standard Serological Pipette	1/100	1mL	Yes	Individually paper-plastic wrapped	200 Pieces/Box, 4 Boxes/Case
J2401	1 mL Standard Serological Pipette	1/100	1mL	Yes	Individually all-plastic wrapped	50 Pieces/Bag, 20 Bags/Case
J2601	1 mL Standard Serological Pipette	1/100	1mL	Yes	Individually all-plastic wrapped	200 Pieces/Box, 4 Boxes/Case

2 mL Standard Serological Pipette




Stretch

Item	Product Name	Graduation/mL	Volume	Sterile	Packaging	Packaging Specification
J2002	2mL Standard Serological Pipette	1/50	2mL	Yes	Individually paper-plastic wrapped	50 Pieces/Bag, 10 Bags/Case
J2202	2mL Standard Serological Pipette	1/50	2mL	Yes	Individually paper-plastic wrapped	200 Pieces/Box, 4 Boxes/Case
J2402	2mL Standard Serological Pipette	1/50	2mL	Yes	Individually all-plastic wrapped	50 Pieces/Bag, 10 Bags/Case
J2602	2mL Standard Serological Pipette	1/50	2mL	Yes	Individually all-plastic wrapped	200 Pieces/Box, 4 Boxes/Case

DATA SHEET


5 mL Standard Serological Pipette



Stretch

Item	Product Name	Graduation/mL	Volume	Sterile	Packaging	Packaging Specification
J2005	5mL Standard Serological Pipette	1/10	5mL	Yes	Individually paper-plastic wrapped	50 Pieces/Bag, 4 Bags/Case
J2205	5mL Standard Serological Pipette	1/10	5mL	Yes	Individually paper-plastic wrapped	100 Pieces/Box, 4 Boxes/Case
J2405	5mL Standard Serological Pipette	1/10	5mL	Yes	Individually all-plastic wrapped	50 Pieces/Bag, 4 Bags/Case
J2605	5mL Standard Serological Pipette	1/10	5mL	Yes	Individually all-plastic wrapped	100 Pieces/Box, 4 Boxes/Case


10 mL Standard Serological Pipette



Stretch

Item	Product Name	Graduation/mL	Volume	Process	Sterile	Packaging	Packaging Specification
J2010	10mL Standard Serological Pipette	1/10	10mL	Stretch type	Yes	Individually paper-plastic wrapped	50 Pieces/Bag, 4 Bags/Case
J2210	10mL Standard Serological Pipette	1/10	10mL	Stretch type	Yes	Individually paper-plastic wrapped	100 Pieces/Box, 4 Boxes/Case
J2410	10mL Standard Serological Pipette	1/10	10mL	Stretch type	Yes	Individually all-plastic wrapped	50 Pieces/Bag, 4 Bags/Case
J2610	10mL Standard Serological Pipette	1/10	10mL	Stretch type	Yes	Individually all-plastic wrapped	100 Pieces/Box, 4 Boxes/Case

DATA SHEET



Unit: mm

Welding

Item	Product Name	Graduation/mL	Volume	Process	Sterile	Packaging	Packaging Specification
J2110	10mL Standard Serological Pipette	1/10	10mL	Welding type	Yes	Individually paper-plastic wrapped	50 Pieces/Bag, 4 Bags/Case
J2310	10mL Standard Serological Pipette	1/10	10mL	Welding type	Yes	Individually paper-plastic wrapped	100 Pieces/Box, 4 Boxes/Case
J2510	10mL Standard Serological Pipette	1/10	10mL	Welding type	Yes	Individually all-plastic wrapped	50 Pieces/Bag, 4 Bags/Case
J2710	10mL Standard Serological Pipette	1/10	10mL	Welding type	Yes	Individually all-plastic wrapped	100 Pieces/Box, 4 Boxes/Case

25 mL Standard Serological Pipette




Unit: mm

Welding

Item	Product Name	Graduation/mL	Volume	Sterile	Packaging	Packaging Specification
J2125	25mL Standard Serological Pipette	2/10	25mL	Yes	Individually paper-plastic wrapped	25 Pieces/Bag, 8 Bags/Case
J2325	25mL Standard Serological Pipette	2/10	25mL	Yes	Individually paper-plastic wrapped	50 Pieces/Box, 4 Boxes/Case
J2525	25mL Standard Serological Pipette	2/10	25mL	Yes	Individually all-plastic wrapped	25 Pieces/Bag, 8 Bags/Case
J2725	25mL Standard Serological Pipette	2/10	25mL	Yes	Individually all-plastic wrapped	50 Pieces/Box, 4 Boxes/Case

DATA SHEET

50 mL Standard Serological Pipette




Unit: mm

Welding

Item	Product Name	Graduation/mL	Volume	Sterile	Packaging	Packaging Specification
J2150	50mL Standard Serological Pipette	1/2	50mL	Yes	Individually paper-plastic wrapped	25 Pieces/Bag, 4 Bags/Case
J2350	50mL Standard Serological Pipette	1/2	50mL	Yes	Individually paper-plastic wrapped	50 Pieces/Box, 4 Boxes/Case
J2550	50mL Standard Serological Pipette	1/2	50mL	Yes	Individually all-plastic wrapped	25 Pieces/Bag, 4 Bags/Case
J2750	50mL Standard Serological Pipette	1/2	50mL	Yes	Individually all-plastic wrapped	50 Pieces/Box, 4 Boxes/Case

100 mL Standard Serological Pipette



Unit: mm

Welding

Item	Product Name	Graduation/mL	Volume	Sterile	Packaging	Packaging Specification
J2100	100mL Standard Serological Pipette	1	100mL	Yes	Individually paper-plastic wrapped	25 Pieces/Bag, 4 Bags/Case
J2300	100mL Standard Serological Pipette	1	100mL	Yes	Individually paper-plastic wrapped	50 Pieces/Box, 4 Boxes/Case
J2500	100mL Standard Serological Pipette	1	100mL	Yes	Individually all-plastic wrapped	25 Pieces/Bag, 4 Bags/Case
J2700	100mL Standard Serological Pipette	1	100mL	Yes	Individually all-plastic wrapped	50 Pieces/Box, 4 Boxes/Case